

Eleviview™

LIFT THAT LASTS

>45
MINUTES
LIFTING



**LONG LASTING
PERFORMANCE**
EXCELLENT
VISIBILITY

Elevi view is a trademark of Cosmo Technologies Ltd.
Distributed by **FUJIFILM Europe B.V.**

SUBMUCOSAL INJECTABLE COMPOSITION – DESIGNED FOR EASY AND SAFE RESECTION PROCEDURES

INTENDED USE

EleviView™ is a submucosal injectable composition intended to use for resection procedures like EMR for a submucosal lift of lesions prior to excision with a snare or endoscopic device.



5 ampoules / pouch (10 ml / ampoule)

PERFORMANCE



- Long-lasting cushion with > 45 min lifting time
- Appropriate for challenging polyps regardless of size, location, or type
- Improved margin visualisation for less potential risk of perforation

EFFICIENCY



- Requires less volume to create cushions
- Fewer reinjections and piecemeal excisions compared to saline that can increase your time saving

VISIBILITY



- Dyed with Methylene Blue for improved visualisation of target lesion margins

SAFE AND EASY



- Pre-mixed, sterile and ready to use
- Low-viscosity emulsion fits through a 23 G needle and is easy to inject

For more information please contact endoscopy_eu@fujifilm.com or