

**INCLUDING
Latex free
tubes &
balloons**

2017/2018

Guidebook DBE Accessories






FUJIFILM
Value from Innovation

DBE Accessories / FTE Accessori es






Article No.	Quantity	Description	EN-450 P5 20	EN- 450T5	EC- 450BI5	EN- 580T	EN- 580XP	EI- 580BT
-------------	----------	-------------	-----------------	--------------	---------------	-------------	--------------	--------------

Attachment tools





ST-01B		1	BS-2 rubber band attachment tool 41A10773020	●	●	●	●	●
ST-05B (2)		20 pcs/pack	Ballon setting tool	●	●	●	●	●
ST-10		1	Rubber band applicator	●	●	●	●	●

















Article No.	Quantity	Description	EN-450 P5 20	EN- 450T5	EC- 450BI5	EN- 580T	EN- 580XP	EI- 580BT
-------------	----------	-------------	-----------------	--------------	---------------	-------------	--------------	--------------

Overtube

TS-1114B (2)			1	Silicone overtube for EN-580XP, sterile, single-use, with expiration date (Contains silicone rubber)				●	
TS-1214B (2)			1	Silicone overtube for EN-450P5, sterile, single-use, with expiration date (Contains silicone rubber)	●				
TS-1314B (2)			1	Silicone overtube for EN-580T/EN-450T5, sterile, single-use, with expiration date (Contains silicone rubber)		●		●	
TS-12140 (2)			1	Overtube for EN-450P5, sterile, single-use, with expiration date (Contains natural rubber latex)	●				
TS-13140 (2)			1	Overtube for EN-450T5 and EN-580T, sterile, single-use, with expiration date (Contains natural rubber latex)		●		●	
TS-13101 (2)			1	Overtube for EC-450BI5 and EI-580BT, sterile, single-use, with expiration date (Contains natural rubber latex)			●		●
BS-4 (2)			10 sets/pack	Endoscope balloon Ø 35 mm, single-use, with expiration date (Contains silicone rubber) (10 pcs balloon + 20 pcs rubber hand) *ST-10 is needed to attach	●	●	●	●	●
BS-2 (2)			10 sets/pack	Endoscope balloon Ø 35 mm, single-use, with expiration date (Contains natural rubber latex) (10 pcs balloon + 20 pcs rubber band)	●	●	●	●	●
DH-14EN (2)			10 pcs/pack	Distal end hood for EN-450P5-20	●				
DH-17EN (2)			10 pcs/pack	Distal end hood for EN-450T5, EC-450BI5, EN-580T and EI-580BT		●	●	●	●
TY-400 PB-20 PB-30			1	Connection tube kit for silicone overtube, PB-20/30 and 450 Series, exchange once every month or once every 10 cases	●	●	●		
TY-500 PB-20 PB-30			1	Connection tube kit for silicone overtube, PB-20/30 and 500 Series, exchange once every month or once every 10 cases			●	●	
TY-04 PB-20 PB-30			1 set/pack	Connecting tube and filter set (clear tube, filter for endoscope, white tube and filter for overtube), exchange once every month or once every 10 cases	●	●	●		
TY-06 PB-20 PB-30			1 set/pack	One-touch-connector set (2 tubes) for EN-580T, EN-580XP and EI-580BT			●	●	●





















DBE Accessories / FTE Accessori es

Article No.	Quantity	Description	EN-450 P5 20	EN-450T5	EC-450B15	EN-580T	EN-580XP	EI-580BT
KDB-004 	10 sets/pack	Double brush and valve cleaning brush, single-use set, channel max. 2.2 mm	●					●
DBP-250 	10 pcs/pack	Double cleaning brush, single-use, channel max. 2.2 mm, length 250 mm	●					●
DBA-250 	10 pcs/pack	Double cleaning brush, single-use, channel min. 2.8 mm, length 250 mm		●	●	●		●
ECP-10040 	10 pcs/pack	Valve brush, single-use for all endoscopes	●	●	●	●	●	●




















Article No.	Quantity	Description	EN-450 P5 20	EN-450T5	EC-450B15	EN-580T	EN-580XP	EI-580BT*
F1FHEW2520280K-S 	1	FTE-Retrieval basket for single-use only, sterile  , hexagonal Ø 20 mm, 4 wires, teflon tube Ø 2.5 mm, length 280 cm, rod, double LL, intern able for guide wire for 0.035 Inch, vice handle				●		
F1FHEW2530280K-S 	1	FTE-Retrieval basket for single-use only, sterile  , hexagonal Ø 30 mm, 4 wires, teflon tube Ø 2.5 mm, length 280 cm, rod, double LL, intern able for guide wire for 0.035 Inch, vice handle				●		
F1LSEW2618290-LL-S 	1	FTE-Lithotripter spiral for single-use only, sterile  , spiral Ø 2.6 mm, length 290 cm, LL-connector				●		
F1SHEW2520280AK-Z-S 	1	FTE-Retrieval basket (Lithotripter basket) for single-use only, sterile  , hexagonal Ø 20 mm, 4 wires, teflon tube Ø 2.5 mm, length 280 cm, rod able to unscrew, double LL, vice handle				●		
F1SHEW2530280AK-Z-S 	1	FTE-Retrieval basket (Lithotripter basket) for single-use only, sterile  , hexagonal Ø 30 mm, 4 wires, teflon tube Ø 2.5 mm, length 280 cm, rod able to unscrew, double LL, vice handle				●		
F1VSEW1210330-S 	1	FTE- Extension rope for lithotripter basket for single use only, sterile  , M1,2 x 330 cm, for assembly of lithotripter coil				●		
F1NSEW1825250Z-S 	1	FTE-Nitinol Polycatcher for single-use only, sterile  , spiral Ø 25 mm, 5 wires, teflon tube Ø 1.8 mm, length 250 cm, with ergonomic handle	●	●	●	●	●	
F1NSEW2330250-S 	1	FTE-Nitinol Polycatcher for single-use only, sterile  , spiral Ø 25 mm, 4 wires, teflon tube Ø 2.3 mm, length 250 cm, with ergonomic handle		●	●	●		

* With a length of 155 cm the EI-580BT provides compatibility with most standard ERCP devices.













DBE Accessories / FTE Accessori es











Article No.	Quantity	Description	EN-450 P5 20	EN-450T5	EC-450B15	EN-580T	EN-580XP	EI-580BT
F2EZTV1805250HP-S 	5 pcs/pack	FTE-Sclerotherapy-needle-set for single-use only, 5 needles, sterile  , needle Ø 0.5 mm, needle length 5 mm, tube Ø 1.8 mm, length 250 cm, with click-needle-in-click-needle-out-mechanism handle, resistant against Histoacryl	●	●	●	●	●	●
F2ZEFJ1806250HP-S 	5 pcs/pack	FTE-Sclerotherapy-needle-set for single-use only, 5 needles, sterile  , inside tube PEEK, needle Ø 0.6 mm, needle length 5 mm, teflon tube Ø 1.8 mm, length 250 cm, with click-needle-in-click-needle-out-mechanism handle, resistant against Histoacryl	●	●	●	●	●	●
F2EZTV2305250-S 	1	FTE-Sclerotherapy needle for single-use only, sterile  , needle Ø 0.5 mm, needle length 5 mm, teflon tube Ø 2.3 mm, length 250 cm, with click-needle-in-click-needle-out-mechanism handle		●	●	●		●
F3CTEW1810250-S 	1	FTE-ERCP-Catheter for single-use only, sterile  , short-filiform, teflon tube Ø 1.8 mm, with mandrin, double LL, length 250 cm, able for guide wire for 0.035 Inch	●	●	●	●	●	●
F3CPEW1810250-S 	1	FTE-ERCP-Catheter for single-use only, sterile  , metal cannula, teflon tube Ø 1.8 mm, with mandrin, DLL, length 250 cm, able for guide wire for 0.035 Inch	●	●	●	●	●	●
F3QBEW1820250-S 	1	FTE-B2-Papillotome for single-use only, sterile  , type loop, teflon tube Ø 1.8 mm, 20 mm cut, double LL, length 250 cm, ergonomic handle with HF-connector 3 mm	●	●	●	●	●	●
F3LQPK0850650X-S 	1	FTE-Wildcat, sterile  , nitinol exchange guide wire antikink, teflon coat striped. Ø 0.035 Inch (0.89 mm), high flexible tip 0.027 Inch (0.69 mm), length 650 cm	●	●	●	●	●	●
F3LQPK0850650XJ-S 	1	FTE-Wildcat, sterile  , nitinol exchange guide wire antikink, teflon coat striped. Ø 0.035 Inch (0.89 mm), high flexible tip 0.027 Inch (0.69 mm), length 650 cm, J-tip	●	●	●	●	●	●
F3QYEW1802250Z-S 	1	FTE-Pre-Cut-Papillotome for single-use only, sterile  , type knife, teflon tube Ø 1.8 mm, knife 0.8 mm x 0.2 mm, knife length 5 mm, double LL, length 250 cm, ergonomic handle with HF-connector 3 mm	●	●	●	●	●	●
F3QTEW1830250-S 	1	FTE-Papillotome for single-use only, sterile  , tip filiform 5 mm, cut 30 mm, teflon tube Ø 1.8 mm, double LL, length 250 cm, ergonomic handle with HF-connector 3 mm	●	●	●	●	●	●

DBE Accessories / FTE Accessori es

Article No.	Quantity	Description	EN-450 P5 20	EN-450T5	EC-450BI5	EN-580T	EN-580XP	EI-580BT
F3YTEW2230250-FR7-5-S 	1	FTE-Y-3 Lumen-Papillotome for single-use only, sterile  , tip filiform, cut 30 mm, teflon tube 2 Lumen Ø 2.3 mm tapered to 1.8 mm, short tapering 11 o'clock model, length 250 cm, handle with HF-connector 3 mm, able for guide wire for 0.035 Inch		●	●	●		●
F3YBEW2225250-FR7-5-15-S 	1	FTE-B-2-Papillotome for single-use only, sterile type loop retrograd, tube 2.3 mm distal tapered to 1.8 mm, 25 mm cut, tip distal 15 mm, DLL, length 250 cm, handle with HF-connector 3 mm, able for guide wire for 0.035 Inch		●	●	●		●
F3YZEW2202250Z-S 	1	FTE-Y-3 Lumen-Needlepapillotome for single-use only, sterile  , teflon tube Ø 2.3 mm tapered tip, needle Ø 0.2 mm, needle length 5 mm, double LL, length 250 cm, handle with HF-connector 3 mm, able for guide wire for 0.035 Inch		●	●	●		●
F4EHPK1825250P-S 	5 pcs/pack	FTE-Polyectomy loop-set for single-use only, 5 loops, sterile  , hexagonal Ø 25 mm, teflon tube Ø 1.8 mm, length 250 cm, handle with HF-connector 3 mm	●	●	●	●	●	●
F4EOPK1810250MP-S 	5 pcs/pack	FTE-Polyectomy loop-set for single-use only, 5 loops, sterile  , oval Ø 10 mm monofil, teflon tube Ø 1.8 mm, length 250 cm, handle with HF-connector 3 mm	●	●	●	●	●	●
F4EOPK1820250MP-S 	5 pcs/pack	FTE-Polyectomy loop-set for single-use only, 5 loops, sterile  , oval Ø 20 mm monofil, teflon tube Ø 1.8 mm, length 250 cm, handle with HF-connector 3 mm	●	●	●	●	●	●
F4EOPK2310250P-S 	5 pcs/pack	FTE-Polyectomy loop-set for single-use only, 5 loops, sterile  , oval Ø 10 mm, teflon tube Ø 2.3 mm, length 250 cm, handle with HF-connector 3 mm		●	●	●		●
F4EOPK2320250P-S 	5 pcs/pack	FTE-Polyectomy loop-set for single-use only, 5 loops, sterile  , oval Ø 20 mm, teflon tube Ø 2.3 mm, length 250 cm, handle with HF-connector 3 mm		●	●	●		●
F4EOPK2330250P-S 	5 pcs/pack	FTE-Polyectomy loop-set for single-use only, 5 loops, sterile  , oval Ø 30 mm, teflon tube Ø 2.3 mm, length 250 cm, handle with HF-connector 3 mm		●	●	●		●
F4EOPK2330250MP-S 	5 pcs/pack	FTE-Polyectomy loop-set for single-use only, 5 loops, sterile  , oval Ø 30 mm monofil, teflon tube Ø 2.3 mm, length 250 cm, handle with HF-connector 3 mm		●	●	●		●

DBE Accessories / FTE Accessori es

Article No.	Quantity	Description	EN-450 P5 20	EN-450T5	EC-450BI5	EN-580T	EN-580XP	EI-580BT
Accessories								
F7GZEW1835250C-S 	1	 FTE-Tripod polyp retriever for single-use only, sterile  , with 3 hooks, teflon tube Ø 1.8 mm, length 250 cm, grasps sink, with ergonomic handle	●	●	●	●	●	●
F7GREW2340250C-S 	1	 FTE-Tripod polyp retriever for single-use only, sterile  , with 3 round hooks, teflon tube Ø 2.3 mm, length 250 cm, grasps sink, with ergonomic handle		●	●	●		●
F7FMEW1820250-S 	1	 FTE-Foreign body forceps for single-use only, sterile  , teflon tube Ø 1.8 mm, flat rat-tooth jaw 1:2 teeth, length 250 cm, grasps sink, with ergonomic handle, capsule retriever	●	●	●	●	●	●
FPU7270 	1	 FTE-Pusher for endoprosthesis for single-use only, sterile  , 7 fr, length 270 cm, with goldring, X-ray		●	●	●		●

Article No.	Quantity	Description	EN-450 P5 20	EN-450T5	EC-450BI5	EN-580T	EN-580XP	EI-580BT
Reusable								
F3CTPK1810250M 	1	 FTE-ERCP catheter, reusable, Ø tube 1.8 mm, length 250 cm, short-fillform	●	●	●	●	●	●
F5RDPK1630260F 	1	 FTE-Double cleaning brush with metal ball, length 260 cm, bristle 8 mm, Ø bristle 4 mm	●				●	
F5RDPK1950260F 	1	 FTE-Double cleaning brush with metal ball, length 260 cm, bristle 10 mm, Ø bristle 5 mm		●	●	●		●
F6BBPK1802250 	1	 FTE-Biopsy forceps, length 250 cm, Ø jaw 1.8 mm, spoon short-round	●	●	●	●	●	●
F6BOPK2304250 	1	 FTE-Biopsy forceps, length 250 cm, Ø jaw 2.3 mm, length jaw 4 mm, Spoon: oval		●	●	●		●

 = Disposable

 = Autoclave



FUJIFILM

FUJIFILM Europe GmbH

Heesenstr. 31, 40549 Düsseldorf, Germany
Tel.: +49 211-50 89 0, Fax: +49 211-50 89 8700
www.fujifilm.eu, endoscopy@fujifilm.eu

Specifications are subject to change without notice. The name FUJIFILM and the FUJIFILM logo are trademarks of FUJIFILM Corporation. All other trademarks shown are trademarks of their respective owners. All rights reserved.
11/18