

IHE Integration Statement

Vendor	Product Name	Version	Date
FUJIFILM Corporation	Fujifilm FDR-2000AWS/FDR-3000AWS/CR-IR363AWS	V6.0	02/04/2015
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Mammography Image Profile(MAMMO)	Acquisition Modality	N/A	
Mammography Acquisition Workflow Profile(MAWF)	Acquisition Modality	N/A	
Patient Information Reconciliation Profile(PIR)	Acquisition Modality	N/A	
Scheduled Workflow Profile(SWF)	Acquisition Modality	Broad Worklist Query Modality Group Case	
Portable Data for Imaging Profile(PDI)	Portable Media Creator	DVD media	
Internet address for vendor's IHE information: http://www.fujimed.com/sub/product_bulletins.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.fujimed.com/sub/product_bulletins.html		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	