

IHE Integration Statement

Vendor	Product Name	Version	Date
FUJIFILM Medical Systems USA, Inc.	Synapse	3.2.1 SR-362	7/15/2011
The product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Manager Image Display Performed Procedure Step Manager		
Patient Information Reconciliation	Image Manager Performed Procedure Step Manager		
Access to Radiology Information	Image Manager Report Reader		
Simple Image and Numeric Report	Report Reader		
Consistent Time	Time Client		
Teaching File and Clinical Trial Export	Export Selector		
Consistent Presentation of Images	Image Manager		
Key Image Note	Image Manager		
Mammography Image	Image Manager		
Portable Data for Imaging	Portable Media Creator		
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer	Basic Patient Privacy Enforcement Basic Patient Privacy Proof	
	Document Source Integrated Document Source/Repository	Basic Patient Privacy Enforcement	
	Document Repository		
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer Imaging Document Source		
Audit Trail and Node Authentication	Secure Application		
Patient Identifier Cross-referencing	Patient Identifier Cross-reference Consumer		
Patient Demographics Query	Patient Demographics Consumer		
Basic Patient Privacy Consents	Content Consumer	Basic Patient Privacy Acknowledgement View	
Internet address for vendor's IHE information: http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/			
Links to Standards Conformance Statements of the Implementation			
DICOM	http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/		
HL7	http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org	