


FUJIFILM

MEDICAL X-RAY FILM

SUPER HR-T 30



*Easy-to-recycle box
for a clean environment*



MEDICAL X-RAY FILM

SUPER HR-T30

A New High-contrast High-resolution Film for General Radiography


Fujifilm's Super HR-T30 is a high-contrast, orthochromatic film that provides high image quality in even the most demanding diagnostic conditions. Super Thin Σ (sigma)-LIC* 30-Grain Technology delivers ultra-high-resolution, fine-grain images under super-rapid (45-second), rapid (90-second), and longer processing cycles.

Offering high-contrast, low noise images that are ideal for general radiography applications, Super HR-T30's enhanced processing stability maximizes diagnostic image quality regardless of processing variability.


This dramatic improvement in stability also makes quality control easier and more efficient even on a daily basis.

* LIC = Localized Image Center


Characteristic Curve



MTF Curve



RMS Granularity Curve



FUJIFILM

FUJIFILM Corporation

26-30, NISHIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8620, JAPAN